



safe&care



UMA PREOCUPAÇÃO AMBIENTAL POR UM FUTURO MAIS VERDE

São cada vez mais evidentes os impactos negativos que o homem provoca na natureza. Agora mais do que nunca, neste período de pandemia do novo coronavírus, é importante criar práticas sustentáveis.

safe&care

MASK

REUSABLE SOCIAL USE 15 WASHES LEVEL 3

Reusable mask designed to promote group protection. These masks are intended to be worn by people when they leave their homes for travel to their workplace, school or business premises, to help protect the community.

Name of manufacturer: VIEIRA POCARGIL SA

Address: Zona Industrial de S.Martinho do Campo, 4830-102 Povia Lanhoso

Trademark or code/Ref.º: PMSLN3 #SAFE&CARE

Lot N.º. | Month | Year of manufacture: LOT1 /MAIO 2020

Raw material composition: 100% POLIPROPILENO

Size identification (child/adult): Adult



FRONT



VERSE

safe&care



15 LAVAGENS USO GERAL (NÍVEL 3)



CONFECTION COMPLIANCE NORM

TEST/NORM: CONFECTION CONFORMITY / MII

RESULTS:

REF mask. PMSLN3 - NIVEL 3

The mask shows conformity of confection.



PRE-TREATMENT NORM

TEST/NORM: HOUSEHOLD WASHING AND DRYING / ISO 6330:2012

RESULTS:

Washing machine used: Type A

Washing program: 6N (60°C)

Total load (kg): 2

Detergent used: ECE + Sodium Perborate +TAED

Drying process: Suspension

Number of washing and drying cycles: 20



AIR PERMEABILITY NORM AND REQUIREMENTS

TEST/NORM: AIR PERMEABILITY / EN ISO 9237:1995

RESULTS:	A	B
median value (l/(m ² .s) ou mm/s):	189	198
average air flow (l/min):	23	24

MINIMUM REQUIREMENTS:

- Greater than or equal to 8 l/min

NOTE: The minimum requirements values are based on the characteristics of surgical masks certified by EN 14683: 2019 TYPE I

A - after 10 cycles of washing and drying

B - after 15 cycles of washing and drying



PARTICLE RETENTION TEST NORM

TEST/ NORM: PARTICLE RETENTION TEST / MI 142/00

RESULTS:

PRC (greater than or equal to 3 µm) (%) - 97

PRC (0,5 µm to 0,7 µm) (%) - 50

Masks level 2, type I (surgical):

PRC (greater than or equal to 3 µm) - Greater than or equal to 95%

PRC (0,5 µm to 0,7 µm) - Greater than or equal to 35%

Level 2 masks, for professionals in contact with the public:

PRC (greater than or equal to 3 µm) - Greater than or equal to 90%

Level 3 masks, for general population:

PRC (greater than or equal to 3 µm) - Greater than or equal to 70%

NOTE: The minimum requirements values are based on the characteristics of surgical masks certified by EN 14683: 2019 TYPE I

B - after 15 cycles of washing and drying



CONSERVATION AND CLEANING

Maximum number of wash cycles in which mask protection performance is guaranteed: 15

Masks should be washed using a normal, complete wash cycle at 60°C with detergent after each use.

The user will need to ensure that the mask is properly dried before use.



INSTRUÇÕES DE USO USER INSTRUCTIONS

This mask should be used as a complement to the recommended protection and hygiene measures and the rules of social distancing, which are fundamental for the control of COVID-19.

For more information, see <https://covid19.min-saude.pt>

The reusable mask should be used correctly.

The following are recommended instructions for the correct use of the mask, in accordance with the DGS Guideline no. 019/2020 of 03/04/2020:

- a) Hygienise your hands with soap and water or an alcohol-based solution before putting the mask on;
- b) Place the mask with the inner face facing the face, and the outer face facing out;
- c) Fit the rigid end of the mask to your nose, and cover your mouth, nose and chin with the mask, making sure there are no gaps between the face and mask;
- d) Do not touch the mask while it is in use; if this happens, the hands must be automatically washed;
- e) It is recommended to use the mask for a maximum of 4 to 6 hours per day and to wash it before wearing it again.



safe&care

geral@safeandcare.pt

Rua do Parque Industrial nº366

4830-103 Campos PVL

PORTUGAL

T: (351) 253 639 280

FAX: (351) 253 632 600