



**safe&care**

IN ANY ACTIVITY OR PROFESSION  
SAFETY MUST BE FACED IN A WAY THAT  
SERIOUS AND RESPONSIBLE

**PREVENTING THE PROPAGATION OF COVID-19**

**safe&care**

# SURGICAL GOWN

## DISPOSABLE

The manufacturer must accompany the gown with a brochure containing information for the end user, in particular information to allow appropriate use of the gown.

Name of manufacturer: VIEIRA POCARGIL SA  
Address: Zona Industrial de S.Martinho do Campo, 4830-102 Povoá Lanhoso  
Commercial brand or code/Ref.ª: BCDP1 #SAFE&CARE  
Lot n°. | Month | Year of manufacture: LOT1 /MAIO 2020  
Composition of raw-materials: 50% POLYESTER 50% POLIPROPILENO  
Size identification (child/adult): Adult M



### CONFECTION COMPLIANCE NORM

TEST/NORM: CONFECTION CONFORMITY / MI

#### RESULTS:

REF mask. BCDP1 #SAFE&CARE

The surgical gown shows conformity of confection.





#### **BURST RESISTANCE**

TEST/NORM: BURST RESISTANCE/ EN ISO 13938-1:1999 - HUMID

TEST/NORM: BURST RESISTANCE/ EN ISO 13938-1:1999 - DRY

#### **BREAKING STRENGTH AND ELONGATION**

TEST/NORM: BREAKING STRENGTH AND ELONGATION (STRIP METHOD)

EN 29073-3:1992 - HUMID

TEST/NORM: BREAKING STRENGTH AND ELONGATION (STRIP METHOD)

EN 29073-3:1992 - DRY

#### **WATER PENETRATION RESISTANCE**

TEST/NORM: RESISTANCE TO WATER PENETRATION (HYDROSTATIC PRESSURE) IN

ISO 811:2018

#### **BACTERIAL PENETRATION RESISTANCE**

TEST/NORM: RESISTANCE TO BACTERIAL PENETRATION - MEDICAL DEVICES / ISO

22610:2006



#### **INSTRUÇÕES DE USO**

The surgical gown must be handled correctly.

The following are recommended instructions for the correct use of the surgical gown:

- a)** Hygiene your hands with soap and water or an alcohol-based solution before wearing the surgical gown.
- b)** Put on the surgical gown with the internal side facing the body, and the external side facing out.
- c)** Adjust the surgical gown to the body, taking into account the cover of the upper and lower limbs and the proper adjustment at the ends, closing all the adjustment systems.
- d)** The removal of the surgical gown should be done so as not to touch the body by turning the suit inside out.
- e)** The surgical gown should be placed in a suitable hospital waste container.
- f)** The hands shall be cleaned again at the end of the use and handling of the surgical gown.

#### **CARE AND MAINTENANCE INSTRUCTIONS**

Products awaiting use must be stored in a dry, enclosed place and protected from potential contamination.

In case of defect or damage do not use the product, replace with another.

#### **PRODUCT LIFETIME**

The product is for single use only and should therefore be marked with the warning phrase: "DO NOT REUSE".



**safe&care**

[WWW.SAFEANDCARE.PT](http://WWW.SAFEANDCARE.PT)

[geral@safeandcare.pt](mailto:geral@safeandcare.pt)

Rua do Parque Industrial n°366

4830-103 Campos PVL

PORTUGAL

T: (351) 253 639 280

FAX: (351) 253 632 600